

EXPERIMENT 6: PHOTOELECTRIC EFFECT

Related Topics

External photoelectric effect, work function, absorption, photon energy, anode, cathode

Principle

A potassium photo-cell is illuminated with light of different wavelengths. Planck's quantum of action, or Planck's constant(h), is determined from the photoelectric voltages measured.

Tasks

To determine Planck's quantum of action from the photoelectric voltages measured at different wavelengths.

Theory and Evaluation

Half of the inside of the high-vacuum photo-cell is a metalcoated cathode. The anular anode is opposite the cathode. If a photon of frequency f strikes the cathode, then an electron can be ejected from the metal (external photoelectric effect) if there is sufficient energy.

Some of the electrons thus ejected reach the (unilluminated) anode so that a voltage is set up between anode and cathode, which reaches the limiting value U after a short (charging) time. The electrons can only run counter to the electric field set up by the voltage U if they have the maximum kinetic energy, determined by the light frequency,

$$hf - A = \frac{m}{2}v^2$$

where A = work function from the cathode surface, v = electron velocity, m = rest mass of the electron. Electrons will thus only reach the anode as long as their energy in the electric field is equal to the kinetic energy:

$$eU = \frac{m}{2}v^2$$

with e = electron charge = $1.602 \cdot 10^{-19}$ As

An additional contact potential ϕ occurs because the surfaces of the anode and cathode are different:

$$eU + \phi = \frac{m}{2}v^2$$


If we assume that A and ϕ are independent of the frequency, then a linear relationship exists between the voltage U (to be measured at high impedance) and the light frequency ϕ :

$$U = -\frac{A + \phi}{e} + \frac{h}{e} f$$

If we assume $U = a + bf$ to the values measured in the figure below we obtain:

$$h = 6.7 \pm 0.3 \cdot 10^{-34} \text{ Js}$$

Literature value: $h = 6.62 \cdot 10^{-34} \text{ Js}$.



Equipment

Photocell, Interference filters, Spectral lamp, Multimeters.

Set-up and Procedure

The experimental set-up is as shown in Fig. 1. The interference filters are fitted one after the other to the light entrance of the photocell.

The measuring amplifier is used in the following way

- Electrometer $R_e \geq 10^{13} \Omega$
- Amplification: 10^0
- Time constant: 0
- Voltmeter: DC 2V

The high-impedance input of the measuring amplifier is discharged via the 'zero' button between measurements.



Fig.1: Experimental set-up for determining Planck's quantum of action.